Traverse rail



Full Room Coverage








- A traverse rail system makes it possible to transfer anywhere in the room
- Parallel rails can be mounted to the ceiling, or wall to wall
- A manual version is easy to use with little effort from the caregiver
- An electric version is available so the entire system can be used independently
- Accessing an adjoining room is possible with the optional traverse connection


Lift & Care systems for people with reduced mobility

Traverse rail





Dimensions (mm)



Technical specifications

Safe work	king load			24	IO kg
Material .		Aluminum	painted	, ELVZ steel,	ABS

For free demonstration or further information

© 10-04-2020 • Handi-Move • v1.0