



HANDI-MOVE

HANDI-SLIDES

Moving without the need to lift



- Low friction for easy movement
- A wide range for a variety of transfers
- Easy to use
- Comfortable and safe
- Allows independent use
- Budget friendly




Handi-Move INTERNATIONAL

SPECIALISTS IN INNOVATIVE AND RELIABLE LIFT & CARE SYSTEMS

TECHNICAL SPECIFICATIONS


HANDI-SLIDE



Art.nr: 3600 (standard)
Art.nr: 3602 (compact)

- Easy handling through low friction.
- Moving people safely in their beds.
- Standard 122 x 71 cm
- Compact 72x 70 cm


HANDI-SLIDE PADDED



Art.nr: 3605 (large)
Art.nr: 3606 (small)

- Padded version of the Handi-slide.
- Ideal for lateral transfers.
- Handi-Slide Padded large 190 x 65 cm (3605)
- Handi-Slide Padded small 72 x 70 cm (3606)


HANDI-SLIDE 1 WAY



Art.nr: 3603 (chair)
Art.nr: 3604 (bed)

- Glides in one direction only.
- Prevents slumping in bed or chair.
- For use in chairs 45 x 40 cm
- For use in bed 80 x 70 cm


HANDI-TRANSFERSLIDE



Art.nr: 3615

- Lightweight plastic construction.
- Easy sliding transfer in sitting position.
- Allows independent use
- Solid and durable.
- Length 666 mm
- Width 200 mm
- Thickness 25 mm
- Weight 0,9 kg


HANDI-TURNTABLE



Art.nr: 3613 (duo 380 mm)

- Easy turning transfer from bed to (wheel)chair or toilet.
- Prevents shuffling feet when turning.
- Diameter 38 cm
- Height 1 cm
- Safe working load 127 kg

HANDI-HANDLINGBELT




Art.nr: 3608 (closed)

Art.nr: 3609 (open)

- Padded belt.
- Aid for standing transfers, repositioning in bed and leg transfers.
- Height 21,5 cm
- Width 62 cm (3609)
- Height 16,5 cm
- Contour adjustable from 41 cm to 114 cm (3608)


HANDI-ROTATINGCUSHION



Art.nr: 3614

- Soft padded turntable.
- Used for turning in sitting position on chair or bed.
- Allows independent use.
- Diameter 38 cm
- Height 0,8 cm

ROLBORD



Art.nr: 3616: standard

Art.nr: 3617: flexible

- 160 cm x 48 cm

For free demonstration or further information :