

HANDI-MOVE

FRICTIONLESS SLIDES

For safe patient handling in simple transfers



SUREHANDS®



- Low friction for easy movement
- Used for a wide range of transfers
- Comfortable and safe
- Facilitates safe patient handling
- Aides in reducing back injuries to staff
- Budget friendly




SUREHANDS®

SPECIALISTS IN INNOVATIVE AND RELIABLE LIFT & CARE SYSTEMS

TECHNICAL SPECIFICATIONS



SLIDING SHEET




Art.nr: 3600 (standard)

Art.nr: 3602 (compact)

- Easy handling through low friction
- Moving people safely in their beds
- Standard 48 x 28 Inch
- Compact 28.4 x 27.5 Inch

SLIDING MAT




Art.nr: 3605 (large)

Art.nr: 3606 (small)

- Padded version of the Sliding Sheet
- Ideal for lateral transfers
- Standard 74.8 x 25.6 Inch
- Compact 28.4 x 27.6 Inch

ONE-WAY SHEET




Art.nr: 3603 (chair)

Art.nr: 3604 (bed)

- Glides in one direction only
- Prevents slumping in bed or chair
- For use in chairs 17.71 inch x 15.74 inch (3603)
- For use in bed 31.49 inch x 27.55 inch (3604)


TRANSFER BOARD



Art.nr: 3615

- Lightweight plastic construction
- Easy sliding transfer in sitting position
- Allows independent use
- Solid and durable
- length 26.22 Inch
- Width 7.87 Inch
- Thickness 0.98 Inch
- Weight 1.98 Lbs


ROTATION DISK



Art.nr: 3613 (duo 380 mm)

- Easy turning transfer from bed to (wheel)chair or toilet
- Prevents shuffling feet when turning
- Diameter 15 Inch
- Depth 0.32 Inch
- Safe Working Load 280 Lbs

SAFETY BELT




Art.nr: 3608 (closed)

Art.nr: 3609 (open)

- Padded belt
- Aid for standing transfers, repositioning in bed and leg transfers
- height 8.46 inch (3609)
- Width 24.40 inch (3609)
- height 6.49 inch (3608)
- Contour adjustable from 16.14 inch to 44.88 inch (3608)


ROTATING CUSHION



Art.nr: 3614

- Soft padded turntable
- Used for turning in sitting position on chair or bed
- Allows independent use
- Diameter 15 Inch
- Height 0.3 Inch

ROLBORD



Art.nr: 3616: standard

Art.nr: 3617: flexible

- 63 Inch x 18.9 Inch

Local Representative: