

Model 2781-2782

FREE-STANDING RAILS SYSTEM

Free-standing rail system which sets up with ease.




- No holes needed in walls or ceiling.
- Easy to set up or relocate.
- Serves as temporary or permanent solution.
- Length of the rails can be adapted and side supports can be placed against the wall.
- Height and width adaptable to the available space.
- Compatible with a ceiling motor or a portable motor.



Lift & Care systems for people with reduced mobility

www.handimove.com

FREE-STANDING RAILS SYSTEM



● Technical specifications

GENERAL

Safe load:240 kg
 Rail:reinforced, aluminum white lacquered
 Side supports: aluminum, gray lacquered

DIMENSIONS

A Length:..... max 5000 mm (choose the rail length based on the available area)
B Height:max 2500 mm (can be shortened based on the available height)
C Length foot: 1000 mm
D Height foot:..... 30 mm



For free demonstration or further information

