

Model 1260-1265



Digital scale

Lifting and weighing at the same time



- Can be used in combination with almost all Handi-Move hoists
- Automatic zeroing
- Large and clear screen
- Minimal loss of height
- 2 years warranty
- Accurate weighing (100 gr)
- Class III

Lift & Care systems for people with reduced mobility

www.handimove.com


Digital scale

● 2 Models

Digital scale Class III Art. #1260
Calibrated digital scale Class III Art. #1265

● Dimensions and specifications

Maximal scale loading: 300 kg
Batteries: 6 x 1,5 V type AAA
Maximum height: 240 mm
Maximum width: 160 mm
Maximum depth: 65 mm



For free demonstration or further information

