

Model 1074-1079



Thigh supports

The key to making the body support self-adjust for you



- Handi-Move offers a variety of thigh supports to make it possible for many different people to enjoy the benefits of the SureHands® body support
- They make it possible to distribute the pressure and the body weight uniformly. They also ensure that the body is in balance during lifting and that the back stays nice and straight.
- Thigh supports fit securely under the thigh and are made of a maintenance friendly rust-resistant material
- The patented thigh support fix for extra comfort

Lift & Care systems for people with reduced mobility


www.handimove.com

Thigh supports

● Versions and dimensions



Art. n°	Description	Color	A [mm]	B [mm]
1077S	Standard child	gray	160	130
1077M	Standard medium	gray	220	180
1077L	Standard adult	gray	310	250
1074S	Standard Slipproof child	black	158	130
1074M	Standard Slipproof medium	black	220	180
1074L	Standard Slipproof adult	black	285	250
1075M	Contour Slipproof	black	220	150
1078M	Contour	gray	220	150
1076S	Muscular dystrophy slipproof child	black	165	130
1076M	Muscular dystrophy slipproof adult	black	220	175
1079S	Muscular dystrophy standard child	gray	160	130
1079M	Muscular dystrophy standard adult	gray	220	180



For free demonstration or further information

