Wall lift



A versatile solution for many situations







- Portable lifting device
- Wall or floor mounted
- Easy to remove and re-use
- Perfect for areas where access is restricted
- Accommodates both the body support and conventional slings
- Extendible boom arm
- Rotating range 180°
- An economical solution

Lift & Care systems for people with reduced mobility



Wall lift

Three possibilities for fixation

Mounted on wall plate(s)

Mounted on standing support(s)

Wall plate + extension for a larger range

38 cm

29 cm

- Dimensions and specifications
- Actuator directly on the mains power with pneumatic handcontrol
 Actuator on 24 V DC transformer with emergency power facilities and electric handcontrol



Safe working load	145 kg	- и	
Wall plate	width	140 mm	
	height	687 mm	
Standing support	outside dimensions	80 x 50 mm	
	max. height	3000 mm	
Boom arm	extendible	extendible	
Weight	16 kg	16 kg	
Safety mechanism	built-in		

Lifing range

