

Model 3473-3478

Seat sling pool/bath PVC

Comfortable in the water



- Sufficient cut-outs for water drainage
- Assumes the form of the bath
- Reinforced seating surface and backrest
- Adjustable loop length

Lift & Care systems for people with reduced mobility



www.handimove.com

Seat sling pool/bath PVC

● Dimensions and specifications

A.....	38 / 14.96 (cm / inch)
B.....	58 / 22.83
C.....	65 / 25.59
D.....	38 / 14.96
E.....	101 / 39.76
F.....	60 / 23.62
G.....	41 / 16.14
H.....	60 / 23.62

Material: PVC, polyester, plastic
 Maximum load: 225 kg / 496 Lbs



For free demonstration or further information