# Handi move

## **Pool Lift Chair**

Comfortably in and out of the water



- Easy to use
- Ergonomic sitting position
- Sturdy plastic design
- Safe, comes standard with a waist band
- Fits in any (shower) wheelchair
- For both children and adults


Lift & Care systems for people with reduced mobility




### **Pool Lift Chair**

### Dimensions and specifications

Α	
В	.450 mm
C	500 mm
D	450 mm
F	475 mm





Chair	100% PVC
Bands	100% Polyester
Band guides	100% Hostaform
Safety clip	100% stainless stee
Weight	5,5 kg
Expected lifespan	10 years
Maximum load	182 kg

For free demonstration or further information

© 12-04-2023 • Handi-Move • v1.0