

Handi-slides

Moving without the need to lift




- Low friction for easy movement
- A wide range for a variety of transfers
- Easy to use
- Comfortable and safe
- Allows independent use
- Budget friendly

Lift & Care systems for people with reduced mobility


Handi-slides

● Handi-slide




Art.nr: 3600 (standard)
Art.nr: 3602 (compact)
 Easy handling through low friction.
 Moving people safely in their beds.
 Standard 122 x 71 cm
 Compact 72x 70 cm

● Handi-slide padded




Art.nr: 3605 (large)
Art.nr: 3606 (small)
 Padded version of the Handi-slide.
 Ideal for lateral transfers.
 Handi-Slide Padded large 190 x 65 cm (3605)
 Handi-Slide Padded small 72 x 70 cm (3606)

● Handi-slide 1 way




Art.nr: 3603 (chair)
Art.nr: 3604 (bed)
 Glides in one direction only.
 Prevents slumping in bed or chair.
 For use in chairs 45 x 40 cm
 For use in bed 80 x 70 cm

● Handi-transferslide




Art.nr: 3615
 Lightweight plastic construction.
 Easy sliding transfer in sitting position.
 Allows independent use
 Solid and durable.
 Length 666 mm
 Width 200 mm
 Thickness 25 mm
 Weight 0,9 kg

● Handi-turntable




Art.nr: 3613 (duo 380 mm)
 Easy turning transfer from bed to (wheel)chair or toilet.
 Prevents shuffling feet when turning.
 Diameter 38 cm
 Height 1 cm
 Safe working load 127 kg

● Handi-handlingbelt




Art.nr: 3608 (closed)
Art.nr: 3609 (open)
 Padded belt.
 Aid for standing transfers,
 repositioning in bed and leg transfers.
 Height 21,5 cm
 Width 62 cm (3609)
 Height 16,5 cm
 Contour adjustable from 41 cm
 to 114 cm (3608)

● Handi-rotatingcushion



Art.nr: 3614
 Soft padded turntable.
 Used for turning in sitting position
 on chair or bed.
 Allows independent use.
 Diameter 38 cm
 Height 0,8 cm

● Rolbord



Art.nr: 3616: standard
Art.nr: 3617: flexible
 160 cm x 48 cm

For free demonstration or further information

