Handi-slings



A sling for every situation







- Complete range meeting every need
- Extremely comfortable sitting or lying posture
- Most slings available with or without head support
- Made of net fabric: also appropriate for bathing
- Seat always remains free
- Usable in combination with every Handi-Move lifting system
- Made-to-measure

Lift & Care systems for people with reduced mobility

lings m

Handi-slings

	without head support	with head support	with head support and loops	comfort
STANDAARD	 Universal use Seat remains free from the coccyx Suspension on 4 or 2 points Padded thigh flaps Tubform 	 Universal use Seat remains free from the coccyx Moonshaped head support Suspension on 4 or 2 points Padded thigh flaps Tubform 	 Universal use Seat remains free from the coccyx Extra loops at the head support for support Padded thigh flaps Suspension on 4 or 2 points 	 Safe and comfortable feeling for users Seat remains free from the coccyx Moonshaped head support Padded thigh flaps Suspension on 4 or 2 points
	two-piece	with head support	Tubform without head support	Tubform thorax
TOILET	 Ideal for toilet transfer Easy to apply Suspension on 2 or 4 points Padded leg strap 	Ideal for toilet transfer Seat and lower back remain free Easy to apply Also available as bathing sling (PVC) For those with limited head control Padded thigh flaps Suspension on 4 or 2 points	 Ideal for toilet transfer Seat and lower back remain free Easy to apply Also available as bathing sling (PVC) Padded thigh flaps Suspension on 4 or 2 points 	Perfect for toilet use Seat and lower back remain free Crossed straps for better support Plain fabric Suspension on 4 or 2 points Padded thigh flaps
	1 leg amputation	double leg amputation	scoliosis	supine
SPECIFIEK	 Specially designed for people with a single amputation Seat remains free from the coccyx Padded thigh flaps Suspension on 4 or 2 points Tubform 	 Specially designed for people with a double amputation Suspension on 4 or 2 points With/without toilet cut-out 	 Special design for extra back support 4 strong ribs Suspension on 2 or 4 points Padded thigh flaps 	 Universal use Transfer to horizontal position Double layer solid fabric Suspension on 4 points

For free demonstration or further information

© 26-02-2016 • Handi-Move • v2.1